



# **DIGITAL IMAGE PROCESSING**

**MEDICAL IMAGE ANALYSIS**

**MOHAMMED YOUSEF  
YASSIN ELASFAR  
SARA ABOALYAZEED**


# INTRODUCTION

- Pneumonia is a significant health problem worldwide, often diagnosed through chest X-rays.
- Developing an automated system to detect pneumonia can aid in early diagnosis and treatment, reducing the burden on healthcare professionals.



# DATASET

- CHEST X-RAYS PNEUMONIA" DATASET FROM KAGGLE.
- INCLUDES 5,863 IMAGES CATEGORIZED INTO NORMAL AND PNEUMONIA CASES




# PROBLEM GAP

- High rates of false positives and negatives in existing models
- Existing Methods doesn't rely on Digital Image processing a lot.

# OUR CONTRIBUTION

- Use of advanced pre-processing techniques to improve image quality and model performance.


# WORK FLOW





# IMAGE PRE-PROCESSING

- Histogram Equalization
- Noise Reduction(Gaussian Filter)
- Sharpening using Laplacian Filter
- Resize
- Normalization




# HISTOGRAM EQUALIZATION

- Computes the histogram of pixel intensities.
- Distributes the intensities more uniformly across the histogram.
- Improves visibility of features in both dark and bright regions.



# NOISE REDUCTION


- Applies a Gaussian filter, which is a type of low-pass filter.
- Convolves the image with a Gaussian kernel.
- The degree of smoothing is controlled by the standard deviation (sigma) of the Gaussian function.






# SHARPENING

- Applies a Laplacian filter to highlight regions of rapid intensity change (edges).
- Can be combined with the original image to create a sharpened effect.
- Helps to enhance important features such as the edges of lung structures.



# RESIZING


- Scales all images to a predefined size (e.g., 224x224 pixels).
- Ensures consistency in input data, reducing computational complexity.
- Maintains the aspect ratio to avoid distortion (if applicable).



# NORMALIZATION



- Divides pixel values by the maximum possible value
- Sometimes subtracts the mean and divides by the standard deviation of pixel values.
- Helps in speeding up the convergence of the CNN during training.

Normalized Image




# CNN

- Input Layer: Accepts resized, normalized chest X-ray images.
- Convolutional Layers:
- Multiple layers to extract features from images.
- Use of filters/kernels to detect edges, textures, and patterns.
- Relu activation functions to introduce non-linearity.
- Pooling Layers:
- Typically max-pooling to reduce the spatial dimensions and computational load.
- Helps in making the model invariant to small translations.
- Fully Connected Layers:
- Flatten the output from convolutional layers.
- One or more dense layers to interpret the features and make final classifications.
- Output Layer:
- Softmax activation function for multi-class classification (e.g., normal vs. pneumonia).




# RESULTS

## MODEL ACCURACY




## MODEL LOSS GRAPH





# RESULTS

## CONFUSION MATRIX





**THANK**  
**YOU**